

Nexo Medico AI - Complete Project Execution Guide

Multi-Tenant Hospital Management & Artificial Intelligence Supercomputer Platform

1. Project Overview

Nexo Medico AI is a Next.js 14 multi-tenant hospital platform with 7 user roles (Super Admin, Hospital Admin, Doctor, Nurse, Pharmacist, Lab Tech, Patient) and specialized AI Clinical & System Security Engines.

2. Step-by-Step Instructions: How to Run the Project

Step 1: Open Terminal in Project Directory

Navigate to your project root folder:

```
cd "c:\Users\yashw\OneDrive\Desktop\Major Project\Project"
```

Step 2: Install Node Dependencies

Install all required packages (Next.js, Prisma, React, Tailwind, Lucide Icons):

```
npm install
```

Step 3: Database Setup & Migration

Push the Prisma schema to create/sync your SQLite database (dev.db):

```
npx prisma db push
```

Step 4: Seed Database with Initial Data

Populate the database with pre-configured hospital tenants, doctors, nurses, patients, and support issues:

```
npm run db:seed
```

Step 5: Launch Next.js Development Server

Start the local application server:

```
npm run dev
```

Step 6: Access Application in Web Browser

Open your browser (Chrome/Edge) and visit:

```
http://localhost:3000
```

Nexo Medico AI - Complete Project Execution Guide

Multi-Tenant Hospital Management & Artificial Intelligence Supercomputer Platform

3. Pre-configured Demo Login Credentials

All accounts use default password: password123

Role	Email Address	Target Route	Scope
Super Admin	superadmin@nexomedico.ai	/admin/dashboard	All Platform Tenants
Hospital Admin	admin@metrohospital.com	/admin/dashboard	Metropolitan Hospital
Doctor	dr.smith@metrohospital.com	/doctor/dashboard	Cardiology / OP
Nurse	nurse.sarah@metrohospital.com	/nurse/beds	Inpatient Wards
Pharmacist	pharma.alex@metrohospital.com	/pharmacy/dashboard	Hospital Pharmacy
Lab Tech	lab.tech@metrohospital.com	/laboratory/dashboard	Clinical Lab & OCR
Patient	john.doe@gmail.com	/patient/dashboard	Personal Health Portal

4. How to Open & Inspect Database (Prisma Studio)

To visually browse, edit, filter, or search all 55+ database tables in your web browser:

1. Open a new terminal in the Project folder.
2. Run: `npx prisma studio`
3. Open browser at: `http://localhost:5555`

5. Useful Commands & Maintenance

- Type Verification:

`npx tsc --noEmit`

- Reset & Re-seed Database:

`npx prisma db push --force-reset && npm run db:seed`